


# Quick Reference



Version E.1



# Table of Contents

Glossary of Terms .....	4
Epic & Alaris Pump Overview .....	5
Alaris Pump Infusion .....	6
Starting a Fluid or Medication .....	7
Changing the Volume to be Infused .....	9
Documenting Volume Infused (Infusion Verify).....	10
Starting a Subsequent Bag .....	12
Dose Change – Modified Order.....	13
Flushes/Carrier Fluids/Keep Vein Open .....	14
IV Piggyback (IVPB) or Secondary Infusion .....	15
Documenting Infusion Complete .....	16
Dual Sign Off .....	17
Dual Sign Off with a Rate Change .....	18
Pump Association and Disassociation .....	19
Back Association .....	20
Assuming Care: Infusions Already Associated .....	22
Stopping an Infusion .....	23
Restarting an Infusion .....	24
Starting and Titrating a Medication .....	25
Delay Feature .....	26
Medication Bolus .....	27
Fluid Bolus .....	28
Documenting Bolus Infusion Completed .....	30
Multi-Dose Administrations.....	31
Insulin Calculator.....	32
OB Considerations.....	33
Syringe Module.....	34
Flushing a Syringe .....	35
Out of Scope Infusions .....	36
ED Narrator.....	37
Using Rover.....	38
Problem Solving .....	40
FAQs .....	43
Notes .....	45

# Glossary of Terms

**Associate the pump:** Link a specific order to a specific pump channel. This is exclusively a 1:1 connection, and this connection is required for the pump data to pull into Epic. Typically, you will do this by scanning the pump barcode when administering an infusion.

**Back associate the pump:** Link an order to a pump channel after the pump has already started and is running. For example, emergent situations where you had to manually program the pump, and you go back into Epic later to create the association.

**Disassociate the pump:** If discontinuing IV fluids or a medication and the channel is no longer needed, then you disassociate the pump channel from the electronic medical record (EMR). This stops the communication between EPIC & the pump.

**Infusion Verify:** an activity, found in the I&O's flowsheet, that allows you to verify (document) rate changes, stops, and volume directly from the pump into the EMR.

**Quick Restart for Alaris Pump:** Power the pump down on the channel, using **Channel Off**, and then touch any other button on the main keypad ("push 5 to stay alive") before 5 seconds and the pump will clear, but not shut off completely (it will remain in an idle state). It will allow you to reprogram from a fresh start without shutting down and booting up again. Very helpful with intermittent IVPBs.


**Subsequent Bag:** an additional bag of the exact same fluid or medication from a **continuous** order in Epic that has already had at least one bag hung.

**Validating Data:** reviewing and accepting data (documenting) in the Infusion Verify activity in Epic. This is similar to vitals devices currently communicate with EPIC. Pump data (i.e., titrations, stops/pauses, volumes) all flows into Infusion Verify. For data to be saved to the patient

record, you must review and accept this information via Infusion Verify.

# EPIC & Alaris Pump Overview

- The purpose of Epic and pump integration is to increase patient safety, decrease pump programming errors, and document true intake volumes and stop times.
- The nurse is always responsible in verifying the standard medication safety checks.
- There are departments that are out of scope and will not use interoperability.
- There are medications and fluids that are out of scope. These medications will not be charted via interoperability. (see pg 32)
- Medications delivered on any pump other than the Alaris pump are out of scope.
- Medications that are not to be run with interoperability will have a message on the MAR (shown below). In those cases, EPIC will not prompt for a pump scan. The pump must be manually programmed using Guardrails library.



- Pump interoperability relies on WIFI to send communications.
  - Green light = pump connected to WIFI
  - No green light = pump not connected to WIFI
    - Is it just the pump?
      - Try another pump in the same room
      - Check that “wireless connection” is enabled (options menu)
    - Is it the WIFI?
      - Multiple pumps tried and no WIFI connected

Biomed = device issue

IT = WIFI issue

# Alaris Pump Infusion



## Infusion—General Information

To use the Alaris pump integration, the following must be true:

- Orders must be in Epic & verified by pharmacy
- Ensure patient's height and weight have been documented.
- The pump must be on the channel select or home screen before sending a new order to the pump.
- The pump must be connected to WIFI




- Make sure you choose the correct profile on the pump, based on the patient you are caring for:
  - Adult – CC & MedSurg
  - OB
  - Oncology/OP INF-Adult
  - Peds < 10 kg
  - Peds 10-30 kg
  - zPeds >30 kg



# Starting a Fluid or Medication

## Starting a Fluid or Medication

1. Scan the patient's armband.
2. Scan the medication barcode.
3. When the Scan the Infusion Pump window appears, in Epic, scan the barcode on the pump channel



**Pump not scanned**

Scan the pump to continue.


Action:

Pump:

Override Reason:



4. Verify Dose, Rate, & VTBI in MAR window
5. Select Primary under Pump Information



**Pump Information**

Pump: LVP 10938662 (Primary) [Change Pump](#) Start: 02/22 0953

Channel Type:

Volume to Infuse:  
1,000 mL (based on order volume)

6. Click Send Details
7. Epic will send information to the pump channel that was scanned. During this process, a Communication Status message will appear in Epic.




**Pump Communication**

Details sent. Waiting for response...  
Pump has 30 seconds to respond.

**Administration Details**

	MAR	Pump
Dose	18 Units/kg/hr	—
Rate	10.7 mL/hr	—
Concentration	100 Units/mL	—
Medication	heparin (generic) in D5W	—
Weight	55.4 kg	—



**Pump Communication**


Review and start pump now  
Pump has 120 seconds to respond.

**Administration Details**

	MAR	Pump
Medication	0.9% NaCl infusion	—
Rate	50 mL/hr	—
Medication	sodium chloride 0.9%	—


# Starting a Fluid or Medication (cont.)

- Verify that the Rate, VTBI (Volume to be Infused) and all other information on the pump matches what is in the Epic order.(see page 9 for how & when to change the VTBI so infusion does not run dry)
- Press Start on the pump.
- The Communication status message in Epic will show pump started successfully.



Ensure you are ready to press Start on the pump to begin your infusion. Take into consideration the pump's location in the room in relation to the computer and if the tubing is spiked/primed before or after you scan the fluid.

- If Epic times out, press Start on the pump and click Try Again in Epic. Must complete this step to complete documentation.





# Changing the Volume to be Infused


## Changing the Volume to be Infused (VTBI)—(after infusion started)

Ensure that you have **pressed Start on the pump first**, to create the connection between Epic and the Alaris Pump.

**DO NOT change the volume to be infused BEFORE starting the pump.**

After starting the infusion, the **Rate** and the **VTBI** can be changed on the pump.

Pump programming is being pulled directly from the order in Epic; the volume to be infused must match the order (ex- 1000 ml of Normal Saline) You must start the pump first, then go back and adjust your VTBI

1. Press **Channel Select**.
2. Press the **VTBI** button. 
3. Using the keypad enter the new volume.
4. Press **Start**.

# Documenting Volume Infused

## Documenting


Using the **Infusion Verify** button on the toolbar OR the **I/O flowsheet** column drop down menu will pull true infused volumes from the pump.

Always check the “Time” field in the Infusion Verify window to ensure you are documenting to the appropriate timeframe.


- Using the Infusion verify button on the toolbar will open the Infusion Verify window at the current time.
  - “Time” field should be changed if you are trying to back chart
- Using the flowsheet column drop down menu will open the Infusion Verify window at the specified column time.

You will always need to use Infusion Verify to pull in the volume from the pump. Do not use the Volume Calculator/Rate Verify All in the flowsheets.

Note: ED can perform this workflow directly from the Narrator using the icon in the graphic below. In the ED Narrator, this is called Pump Rate Verify



1. Access the I/O flowsheet.
2. Click Infusion Verify on the toolbar OR I/O flowsheet column drop-down menu.
3. Ensure the “Time” field contains the desired time you want to document.



# Documenting Volume Infused (cont.)

4. Choose the rates you want to verify (document) by making sure a checkmark is next to the correct rate row.

Infusion Verify

Date: 10/06/2024 Time: 0758 Now

Select All Collapse All

0.9% NaCl infusion Start: 10/06/21 0800 Infused 46.37 mL

Administrations to file: 4  
From 10/06 07:11:37 to 07:58:00

MAR Events			Volumes			
File	Time	Event	Dose	Rate	Interval	Volume
	10/06 07:11:37	✓ New Bag		100 mL/hr	07:11:37 – 07:11:43	0.17 mL
<input checked="" type="checkbox"/>	07:11:43	Verify Only		100 mL/hr	07:11:43 – 07:16:32	8.02 mL
<input checked="" type="checkbox"/>	07:16:32	Paused		0 mL/hr	07:16:32 – 07:35:05	0 mL
<input checked="" type="checkbox"/>	07:35:05	Restarted		100 mL/hr	07:35:05 – 07:55:18	33.68 mL
	07:55:18	Verify Only		100 mL/hr	07:55:18 – 07:58:00	4.5 mL
<input checked="" type="checkbox"/>	07:58:00	Verify Only		100 mL/hr		

✓ Total Volume ≈ 46.37

Accept Cancel

5. Verify the volumes you are filing are accurate for each infusion.

✓ Total Volume ≈ 46.37

6. Click Accept to verify and file into the flowsheet. You will have to scroll to the bottom to activate the Accept button.

If you do not complete this step, you have NOT documented.

7. The data is now filed in the flowsheet and MAR.


# Starting a Subsequent Bag

## Subsequent Bags


When the volume to be infused is complete, the pump beeps and drops to a KVO rate.

To hang a new bag:

1. Scan the patient's armband.
2. Scan the fluid or medication barcode.
3. When the Scan the infusion pump window appears, scan the barcode on the same pump channel as the current bag.



4. When the Medication Administration window opens, click Send Order Details.



You may have to scroll to the bottom of the window to activate the Send Details to Pump button.

5. Verify the Epic order and the pump exactly match.
6. Press Start on the pump.

**Note: When sending order details, the medication remains running. This step does not pause the infusion**

# Dose Change – Modified Order

## Dose Change - Modified Order


If the order is *modified*, there is no need to rescan the medication or pump. Change the dose on the pump and validate the rate change in the flowsheet.

1. Change the Dose on the pump FIRST.
2. Press Start
3. Open the I&O flowsheet.
4. Click the Infusion Verify button on toolbar OR I/O flowsheet column drop down menu.
5. Ensure you are verifying the correct medication or fluid, which is indicated with a check mark in the File Total Volume box.
6. Uncheck any values you do not want to verify at this time.

*(Reasons to uncheck: pauses, patient bends arm, kink in tubing)*

File	Time	Event	Dose	Rate	Interval	Volume
	10/06	07:11:37	✓ New Bag			
<input checked="" type="checkbox"/>		07:11:43	Verify Only	100 mL/hr	07:11:37 – 07:11:43	0.17 mL
<input checked="" type="checkbox"/>		07:16:32	Paused	100 mL/hr	07:11:43 – 07:16:32	8.02 mL
<input checked="" type="checkbox"/>		07:35:05	Restarted	0 mL/hr	07:16:32 – 07:35:05	0 mL
		07:35:05	Restarted	100 mL/hr	07:35:05 – 07:55:18	33.68 mL
		07:55:18	Verify Only	100 mL/hr	07:55:18 – 07:58:00	4.5 mL
<input checked="" type="checkbox"/>		07:58:00	Verify Only	100 mL/hr		

7. Click Accept to file.



**Note:** If you have not had a chance to Infusion Verify (the previous dose) before the order is modified, you may get a MAR warning the next time you Infusion Verify.

# Flushes/Carrier Fluids/Keep Vein Open

## Flush and Carrier Fluids

Approved standing orders will be available for any flushes/carrier fluids/keep vein open infusions. This order will appear as a “range” order set, which gives nursing the autonomy to adjust the rate as needed.

Order Set: “Flush and Carrier Fluids”.

NOTE: If you have multiple carrier fluids/”KVO” infusions running, you will need an order for each one.

- The Flush and Carrier Fluids order can be used for infusion solutions of NS, 5%D, and LR. You choose what you need.
- There are up to three flush/carrier infusion lines for each type of solution
  - The flush/carrier infusion lines are indicated as A, B, and C.
  - This does not correlate with pump channels
  - A is one line, B is two lines, and C is three lines.
- There is an Adult section and Pediatric section, only choose one


The screenshot displays a software interface for configuring an order set. At the top, the word 'FLUSH' is visible. Below it, a section titled 'Order Sets & Panels' shows a list with 'Flush and Carrier Fluids' selected. A pop-up window titled 'Orders' shows a tree view with 'Flush and Carrier Fluids' expanded to show 'Adults' and 'Peds' options. A detailed view of the 'Adults' section shows three checkboxes: 'Adult - Flush Carrier Solution - PRN', 'Adults - Flush Carrier Solutions - Continuous', and 'Flush via IV Push - Adults'. The 'Peds' section is collapsed. At the bottom, a detailed view of the 'Adults' section shows a checked checkbox for 'Adult - Flush Carrier Solution - PRN' with a red exclamation mark icon, and an unchecked checkbox for '0.9% NaCl infusion - carrier fluid A (\$3.00/day)'. Below this, the text 'As needed (PRN), at 1-35 mL/hr, post IV medication adm' is visible.

# IV Piggyback OR Secondary Infusion


## IV Piggyback (IVPB) or Secondary Infusion

**Remember:** A primary bag or flush bag should be started first and Associated to an order. Now you are ready to start a secondary / piggyback infusion.

1. Scan the patient's armband.
2. Scan the barcode on the medication.
3. When the Scan the Infusion Pump window appears, scan the barcode on the pump channel that your primary is associated with.



4. Verify MAR admin instructions
5. If using Mini Bag +- scan fluid component now
6. Choose Piggyback in the Infusion Pump Setup window



7. Click Send Details to Pump button located at the bottom of the MAR Administration window
8. **Important: Press the Secondary button on the pump.**
9. Press Next on the pump.
10. Verify that the Rate, VTBI, and all information on the pump matches the Epic order.
11. Press **Start** on the pump.

# Documenting Infusion Complete

## Infusion Complete

Once the infusion is completed, you MUST validate the volume(s) in the I/O flowsheet to document that the infusion is completed.

1. Access the I/O flowsheet.
2. Click on the Infusion Verify button on the toolbar or I/O flowsheet column drop down menu. The Infusion Verify window will open.
3. Ensure the “Time” field contains the desired time you want to document.
4. Locate the fluid/medication that you want to verify.
5. Notice the zero rate and the Stopped action. Select the magnifying glass to change the action to Paused if needed.

The screenshot shows the 'Infusion Verify' window. At the top, the date is 10/6/2021 and the time is 0905. The infusion is identified as '0.9% NaCl infusion' starting at 10/06/21 0800. The total volume infused is 154.17 mL. A red box highlights the infusion name. Below this, a table lists MAR events and volumes. The 'Stopped' event at 09:02:42 has a rate of 0 mL/hr, which is also highlighted with a red box. The 'Action' dropdown menu is set to 'Stopped' and is also highlighted with a red box. At the bottom, there are buttons for 'Keep Association', 'Disassociate', 'Go to Next', 'Accept', and 'Cancel'. A red warning icon indicates that documentation is missing.

File	Time	Event	Dose	Rate	Interval	Volume
	10/06	07:11:37	✓ New Bag	100 mL/hr	07:11:37 – 07:11:43	0.17 mL
✓	07:11:43	Verify Only	100 mL/hr	100 mL/hr	07:11:43 – 07:16:32	8.02 mL
✓	07:16:32	Paused		0 mL/hr	07:16:32 – 07:35:05	0 mL
✓	07:35:05	Restarted	100 mL/hr	100 mL/hr	07:35:05 – 09:02:42	145.98 mL
✓	09:02:42	Stopped		0 mL/hr	09:02:42 – 09:05:00	0 mL

6. Click Accept.

Note: ED can perform this workflow from the ED Narrator using the icon in the graphic below.







# Dual Sign-off

## Completing a Dual Sign Off


1. Scan the patient's armband.
2. Scan the barcode on the medication.
3. When the Scan the infusion pump window appears, scan the barcode on the pump channel.




4. Select Primary in the Pump Information section and click Accept.



5. Both nurses Verify the medication administration details.
6. Click Send Details.





7. Pump: Verify pump matches order. Click Start
8. EPIC: Complete Dual Sign Off



# Dual Sign Off with Rate Change

## Dual Signoff with a Rate Change


1. Change the rate on the pump.
2. Press Start on the pump.
3. Access the I/O flowsheet & click the Infusion Verify button on the toolbar or I/O flowsheet drop down menu. 
4. Ensure the “Time” field contains the desired time you want to document.
5. Choose the rates you want to verify (document) by making sure a checkmark is next to the correct rate rows.



The screenshot shows the 'Infusion Verify' window. At the top, there are fields for 'Date' (10/06/2024) and 'Time' (0942). Below this, a medication entry is shown: 'heparin 25,000 units in dextrose 5% 250 mL (100 units/mL) infusion HIGH INTENSITY nomogram - OHS' with a start time of 10/06/21 0700 and an infused volume of 33.63 mL. A section titled 'Administrations to file: 3' shows a time range from 10/06 08:06:22 to 09:42:00. Below this is a table with columns for 'File', 'Time', 'Event', 'Dose', 'Rate', 'Interval', and 'Volume'. A red box highlights the 'File' column, which contains checkmarks for several rows. At the bottom right, there is a 'Total Volume' field showing 33.63 mL and buttons for 'Accept' and 'Cancel'.

File	Time	Event	Dose	Rate	Interval	Volume
	10/06 08:06:22	New Bag	18 Units/kg/hr	10.7 mL/hr	08:06:22 - 08:06:31	0.03 mL
<input checked="" type="checkbox"/>	08:06:31	Verify Only	18 Units/kg/hr	10.7 mL/hr	08:06:31 - 09:32:51	15.38 mL
	09:32:51	Bolus Start			09:32:51 - 09:32:51	
	09:32:51	Bolus	250 Units/kg/hr	153.3 mL/hr	09:32:51 - 09:39:51	17.84 mL
	09:39:51	Bolus End			09:39:51 - 09:39:51	
<input checked="" type="checkbox"/>	09:39:51	Verify Only	18 Units/kg/hr	10.7 mL/hr	09:39:51 - 09:40:18	0.08 mL
	09:40:18	Verify Only	10 Units/kg/hr	10.7 mL/hr	09:40:18 - 09:42:00	0.3 mL
<input checked="" type="checkbox"/>	09:42:00	Verify Only	18 Units/kg/hr	10.7 mL/hr		

6. The Dual Signoff window will open; have the second nurse sign off the medication. Click “OK”



A dialog box titled 'Dual Signoff Second User:' with two input fields: 'User ID:' and 'Password:'.

To perform a **Dual Sign Off Bolus** see bottom of page 27


# Pump Association and Disassociation

## Pump Association

This is the linking of a specific order to a specific pump channel. This is exclusively a 1:1 connection, and this connection is required for pump data to pull into Epic. Typically, you will do this by scanning the pump barcode when administering an infusion.

To see all the pump associations:

1. Click the **MAR**.
2. Click the **Pumps** button on the toolbar.



This will show all the fluids/medications that are associated with the channels.


## Disassociations

If discontinuing IV fluids or a medication and the channel is no longer needed, then you disassociate the pump channel from the electronic medical record (EMR). This stops the communication between EPIC & the pump.

**Disassociate does NOT delete the information**

**Disassociate does NOT document a STOP**

Disassociation can be done through Infusion Verify.




Turning the pump “off” does not disassociate.

# Back Association


## Back Association - from an out-of-scope department

When a patient has been transferred from an out-of-scope area such as the OR into an in-scope area such as PACU and the infusion is already running on the pump, you will need to associate the pump.


1. Scan the patient armband.
2. Scan the barcode on the medication.
3. When the product scan window opens, choose the action  
Back Association. (NOT new bag) and press Accept



4. In Pump Information section, click Associate Pump. Then scan the pump.
5. Choose Primary infusion.



Important: Choose the time the you have taken responsibility of the patient on the arriving unit as the Association Start Time. Choose "Use original time" since this is the time that you are taking over patient care.



6. Click Accept

Note: you will NOT follow the "New Bag" process. You are only needing to associate that medication to that specific module.


# Back Association (cont.)

## Back Association – Infusion you initiated

In an emergency or in downtime, you may need to hang fluids/medication on your patient before there is an order in Epic. You can manually program the pump as you currently do, and then Back Associate once the order is in Epic.

### How to Back Associate:

1. Scan the patient armband.
2. Scan the medication barcode.
3. Leave the MAR “action” as New Bag. You are stating that YOU initiated this infusion and hung the first bag.
4. When the Scan the infusion pump channel window appears, scan the barcode on the pump channel.
5. Verify MAR dose & rate
6. Choose Primary in the “Infusion Pump Setup window”.
7. In the Pump Information window document the time you actually started the fluids (this is very important) DO NOT choose “Use original time”



8. Click **Accept**.

You will NOT have to “send details to pump”

# Assuming Care: Infusions Already Associated

You've assumed care of a patient (shift change or admit), and the infusions have already been associated. However, the person who last cared for the patient did not perform an Infusion Verify.

1. Note the time you assumed care of the patient
2. Go into the I/O Flowsheet
3. Access Infusion Verify and ensure the "Time" field contains the desired time you want to document.
4. Uncheck the "big box" first – this unchecks data pulled under this med -helpful if there is a lot of data
5. Ensure the time you are claiming is "checked"

Infusion Verify

Date: 10/6/2021 Time: 0847 Now

Select All Collapse All Got It

When you clear a MAR event check box, the event doesn't file to the MAR, but the volume for that event still counts toward the total volume to file.

0.9% NaCl infusion Start: 10/06/21 0800 Infused 0 mL

Administrations to file: 1  
From 10/06 07:11:37 to 08:47:00

File	Time	Event	Dose	Rate	Interval	Volume
	10/06 07:11:37	New Bag		100 mL/hr	07:11:37 - 07:11:43	0.17 mL
<input type="checkbox"/>	07:11:43	Verify Only		100 mL/hr	07:11:43 - 07:16:32	8.02 mL
<input type="checkbox"/>	07:16:32	Paused		0 mL/hr	07:16:32 - 07:35:05	0 mL
<input type="checkbox"/>	07:35:05	Verify Only		100 mL/hr	07:35:05 - 08:45:18	116.98 mL
<input type="checkbox"/>	08:45:18	Verify Only		100 mL/hr	08:45:18 - 08:47:00	2.83 mL
<input checked="" type="checkbox"/>	08:47:00	Verify Only		100 mL/hr		

Total Volume: 0 Revert to calculated (128 mL)

Accept Cancel

6. Change volume to "0"

7. Click Accept

Remember: Out of scope areas will continue to chart medications and I/O's as they currently do

The next time you Infusion Verify – only the data from your last pump rate verified will pull over- think "Infusion Verify sandwich"


# Stopping an Infusion

## When to Stop an Infusion


- Discontinued Medication (Action: Stopped)
- One-time orders (Action: Stopped)
- Allergic reaction (Action: Stopped)

## Stopping an Infusion

1. Stop the fluid or medication on the pump.
2. Access the I/O flowsheet.
3. Click Infusion Verify button on the toolbar or the I/O flowsheet column drop down menu.



4. Choose the rates you want to verify (document) by making sure a checkmark is next to the correct group.
5. Notice the zero rate and the options “Stopped” and “Paused”. Make sure Stopped is highlighted.



After the infusion is complete, disassociate the pump from that order.

6. Click the Disassociate button, next to the 0 (zero) rate in the Infusion Verify window.
7. Click **Accept**.


# Restarting an Infusion

## Restarting a Medication

(Ex: titratable medication -briefly titrated off – needs to be restarted)

If the brain was powered down, (medication was temporarily stop/disconnected from the patient) you can restore the most recent settings and infusions by re-starting using **Restore** on the pump.

1. Turn on the pump. (If pump is already “on” skip to step 4)
2. Choose No for New Patient.
3. Press Yes to continue with same profile.
4. Press the pump channel button for the fluid that you want to continue.
5. Push Restore on pump.
6. Verify the Epic order matches the pump settings and press Start on the pump.
7. Access the I/O flowsheet.
8. Click Infusion Verify button on the toolbar or the I/O flowsheet drop down menu.
9. Choose the rates you want to verify (document) by making sure a checkmark is next to the correct group. **The action should be Restarted.**




10. Click **Accept**.



# Starting and Titrating a Medication

## Starting a Titratable Medication

1. Scan the patient's armband.
2. Scan the medication barcode.
3. When the Scan the Infusion Pump window appears in Epic scan the barcode on the pump channel.




4. When the Medication Administration window opens, make sure the ordered dose range is correct.



5. Enter the initial ordered dose.



6. In Epic, select Primary in the Pump Information setup window



7. Click Send Details.
8. Verify that the Rate, VTBI and all other information on the pump matches what is in the Epic order.
9. Press Start on the pump.

# Starting and Titrating a Medication (cont.)

## Titrating a Medication (Non-Emergent Situation)


1. Verify titration order parameters in the MAR
2. Adjust dose on the pump as needed per order
3. To document go to I/O flowsheet
4. Click Infusion Verify button on the toolbar or the I/O flowsheet column drop down menu and document actions
5. Click Accept.

(If the MAR Warning window populates, select “Titratable Drug”

This should not happen if the dose is entered as a range order)

## Utilizing the “Delay” Feature

Some departments utilize the “delay” feature on the pump. If there is a reason to delay your pump, ensure you are utilizing the “Delay For” and NOT the “Delay Until”.



# Medication Bolus

## MEDICATION Bolus from an Existing Bag

The Medication must already be infusing on the pump.

1. Scan the patient's armband
2. Scan the medication barcode

The **"Select An Order"** window opens and displays the base order and the bolus.



3. Select the correct bolus order and verify correct dose in MAR. Click Accept.

## On Pump:

1. Press the Channel Select button.
2. Press Bolus.
3. Program the pump per order (dose, VTBI, duration, etc)
4. Press Start.
5. Click Accept in EPIC & complete dual sign-off



When the Bolus is finished, the pump will return to the original rate.

## Dual Sign-off Bolus – from an existing bag

Additional notes:



- Follow above steps in administering the bolus.
- No pump scan required, therefore you will not "send order details to pump".
- Manually program the bolus.
- You will be prompted for the dual sign off.

# Fluid Bolus

## FLUID Bolus from Existing Bag

The fluid must already be infusing on the pump. Consider the volume remaining in the bag. Do you need a separate bag? Or is there enough in the bag that is currently infusing?

1. Scan patient's armband.
2. Scan medication barcode.
3. If given choice to choose order: Select the Bolus order.
4. Click Accept.
5. DO NOT scan the pump\*
6. Change the "Action" to "Bolus from Bag"
7. Click Continue Administration.
8. View MAR to verify bolus order
9. Click Accept in MAR



## On the pump:

10. Press Channel Select.
11. Program your bolus rate and VTBI from the order by adjusting your current primary infusion settings.
12. Press Start on the pump.

When the Bolus is finished, the pump will drop down to KVO. You will need to assess the patient and reprogram the pump for the rate and VTBI from the primary infusion order.

When performing an Infusion Verify, bolus volume is documented under the CONTINUOUS INFUSION NOT the bolus order. Do NOT use the "calculator" under the bolus in the I/O flowsheet to calculate the volume.

lactated ringers infusion			
Start: 11/06/20 1530			
Rate		250	
Volume (mL)			500
lactated ringers bolus 500 ml			
Start: 11/06/20 1530			
Dose			
Volume (mL)			
Rate			

Bolus dose of 500ml documented under continuous LR infusion (via Infusion Verify)

Do NOT use volume calculator under bolus order

# Fluid Bolus (cont.)

## **Fluid Bolus from a New Bag (utilizing a separate channel)**

(not enough volume in current bag – hanging a separate bag on a separate channel to give bolus)

This will be the same workflow as hanging/starting a new primary bag.

See page 7 for “Starting a Fluid or Medication”

## **Fluid Bolus from an Existing Bag (but not enough volume in current bag)**

(Plan to use same pump as your continuous infusion, but the bag needs to be changed in order to give the bolus amount)

See page 12 for “Subsequent Bag”

Then follow page 28 - “Fluid Bolus from Existing Bag”

# Documenting Bolus Infusion Complete

## Documenting Bolus Infusion Completed

Once the bolus is completed, you will need to validate the infusion intake and document the infusion is completed via **Infusion Verify**.

1. Click Infusion Verify button on the MAR or the I/O flowsheet.

When the Bolus is finished, the pump will drop down to KVO. You will need to assess the patient and reprogram the pump for the rate and VTBI from the primary infusion order.

**Volume is documented under the CONTINUOUS INFUSION NOT the bolus order. Do NOT use the "calculator" under the bolus in the I/O flowsheet to calculate the volume.**

lactated ringers infusion			
Start:	11/06/20 1530		
Rate		250	
Volume (mL)			500
lactated ringers bolus 500 ml			
Start:	11/06/20 1530		
Dose			*500 mL
Volume (mL)			
Rate			250

Bolus dose of 500ml documented under continuous LR infusion

Do NOT use volume calculator under bolus order

# Multi-Dose Administrations

## Orders Requiring Multiple Doses

*Example: Provider order 40 mEq Potassium IV and you pull 4 bags of 10 mEq Potassium IV.*

- Remove full dose from the pyxis
- Scan each dose (bag) as you administer them (do NOT scan all IV bags at one administration time)
- A MAR warning will appear, noting that the entered rate/dose does not match the ordered rate/dose
  - Choose Multi-Dose Administration
  - Verify that the dose and rate are correct before sending details to pump.
- When the first dose is completed, repeat these steps to properly document and administer the next dose

**Rate entered is less than rate ordered**

dextrose 5 % and 0.45 % NaCl with KCl 40 mEq infusion  
Entered: 10 mL/hr  
Ordered: 100 mL/hr

Cancel and re-evaluate the rate.

Override Reason:

Title	Number
Bolus from Bag	57
Dose ordered not in med cabinet	22
Imaging Contrast Dose Adjustment	195
Injection stopped due to extravasation	70
<b>Multi-Dose Administration</b>	309
Patient had allergic reaction	69
Patient Refused	11
Patient unable to tolerate	68
Pump in KVO mode	303
Titratable drug	23

# Weight-Based Insulin Infusion Calculator

## Insulin Infusion Calculator-Initiation

1. Scan the patient, medication, and IV pump channel
2. Select **INITIATE** action then enter **the current Glucose** level to generate the calculated dose. Ensure the following are entered in the Med Admin details window located above the calculator” MAR action= **New Bag**, Dose field = **calculated dose value**. Select the pump information/channel type = **Primary**
3. Select the **Order** button on the MAR. Confirm MAR Details are displaying on the infusion pump. You will be prompted to complete a **dual signature** on the MAR.

## Insulin Infusion Calculator Dose Adjustments

1. Go to the MAR and select “**Rate/Dose Change**” to open the Medication Administration Window
2. Use the Insulin Infusion Calculator
  - Select **Calculate Dose** action then enter **current POCT Glucose** to generate the newly calculated dose
3. **GO TO THE PUMP**- Manually change the dose on the pump if the calculated dose is different.
  - If the dose stays the same change MAR action to **No Rate/Dose Change**
4. Click Accept on the MAR and complete dual sign off
5. Perform an infusion verify the next time you document Intake/Output

Because the infusion calculator lives in the MAR, the calculated new dose is done in the MAR administration window. Make the dose change on the pump PRIOR to clicking Accept in Epic (where dual sign is performed).



# OB Considerations

## Magnesium infusion and Bolus

You must start the infusion using interoperability then administer the bolus using the bolus soft key.

**\*\*Don't forget to document the bolus administration on the MAR\*\***

## Oxytocin titrating for a laboring patient

1. Scan the patient, medication, then scan the IV pump channel
2. Validate information in the Epic administration window
3. Click **Send Details**
4. Validate the information on the Alaris pump and select Start

## After placenta delivery

1. Select Channel Select for Oxytocin infusion
2. Select **Bolus** soft key to administer
3. Alaris pump has bolus dose set to 10 u over 30 minutes
4. Confirm your dose
5. Select **Start**
6. Manually document Bolus administration on the MAR Bolus order

## Post-Partum Oxytocin dose (non-titrating)


1. **Quick restart** the pump
2. Choose Guardrail Fluids
  - Select K-O
  - Select oxytocin\_postpartum
  - Manually program the rate and VTBI
  - Press Start
3. Back associate the oxytocin infusion from the MAR using a **"Back Association"** MAR action

# Syringe Module


## Syringe Module

Prepare the medications and tubing as you currently do. Then follow the steps below to continue with the integration into Epic.

1. Scan the patient's armband.
2. Scan the medication barcode.
3. When the Scan the Infusion Pump window appears in Epic scan the barcode on the syringe pump channel.



4. When the "Medication Administration" window opens verify that the order is correct.
5. Click Send Details.
6. Install the syringe into the pump at this time- the pump will tell you when to complete this action.
7. Select the syringe size.
8. Verify the medication details. Click Next.
9. Press Start




# Flushing a Syringe Module

## **Adding a FLUSH syringe after the medication has completed:**

(You must have an order(s) for flushes. They will be charted in the MAR, but the “order details” will NOT be sent to the pump.)

1. Change the syringe to a flush syringe.
2. Press channel select.
3. “Confirm” syringe size.
4. Press Restore on the pump.



**This action will restore the previous medication that was programmed.**

5. Press Next.
6. Change the VTBI to the desired flush volume you wish to instill.
7. Press Start.

**Note: In the I/O flowsheet the flush volume will be included in the medication volume. The MAR will be the source of truth for amount of medication given.**

## **Underfilled Syringe/Line primed with medication:**

- You will have to adjust your VTBI
- You will have to adjust your rate to match the ordered rate
- If EPIC sends a VTBI that is greater than what the pump detects, the populated VTBI will be blank – which requires a manual entry

# Out of Scope Medications

An “out of scope” infusion must be manually programmed on the pump. Meaning, you will need to go into the library, find the infusion, and enter in the details from the order in Epic.

Remember, if you are manually programming the pump you will need to manually document any changes in the MAR and enter your totals into the I&O's Flowsheet.

Below is a list of some common “**Out of Scope**” **Infusions**:

- Blood products
- PCAs
- Epidurals
- Any medication NOT run on an Alaris pump
- Medications administered during a Code Blue
- Penicillin G
- Investigational drugs that are not in the drug library


There will be a notation on the MAR for out-of-scope medications. This is what you will see on the MAR:

NeoFax  
Pediatrics  
Education/Micromedex (MEFS)  
Ask Watson

Order ID: 238986024

disposal code: DEA


Infusion Pump Message: Manually program pump for this medication

Click to see more details 


# ED Narrator

All steps required for interoperability can be completed through the ED Narrator.

- Mini MAR can be accessed on the ED Narrator




- Once an infusion has been started using interoperability, there will be a Pump Rate Verify icon on the mini MAR.




- The Pump Rate Verify window will open when the icon is selected.

NOTE: If all infusions are in the same category, they can be verified together by selecting the icon next to the category name. However, If infusions fall into different categories, each category will need to be verified separately.




NOTE: If using the MAR activity, you will document volumes and actions using Infusion Verify, not Pump Rate Verify



# Using Rover

Rover enabled devices can be used to administer IV fluids and medications using interoperability.

1. Go to the MAR.
2. Create “New Admin”.
3. Scan patient armband.
4. Select the medication or fluid. It then appears in the “cart”.
5. Click Accept.
6. Scan the medication/fluid barcode.
7. Scan the pump when prompted.
8. Select Primary or Piggyback when prompted.  
Then click Next.



Note: You can NOT Infusion Verify via Rover. Infusion Verify must be done using a computer.

## **Back Association from an out-of-scope area (you did not start)**

1. Scan the patient armband and medication barcode.
2. Click “Admin Details” and choose the action Back Association. (NOT new bag)
3. Use the back arrows to access the medication in the Cart
4. Once in the Cart, click “Accept” to scan the pump channel.
5. Select the Channel Type and the Association Start. Choose “Now” since this is the time that you are taking over patient care.

# Using Rover

6. Click Next. The administration is now saved.

## **Back Associate an infusion that you started earlier (during downtime or in an emergency)**

1. Scan the patient armband and medication barcode.
2. Click “Admin Details” and keep action as New Bag.
  - a. You are stating that YOU initiated this infusion and hung the first bag.
3. Select the Dose field. Input the dose currently running on the pump and
4. Use the back arrows to access the medication in the Cart.
5. Once in the Cart, click “Accept” to scan the pump channel.
6. Select the Channel Type and the Association Start. Choose the date and time you started the infusion. Do not choose “Now”.

The screenshot shows a mobile application interface for Rover. At the top, the status bar displays the time 1:53, signal strength, 5G, and battery level. Below the status bar, there are three navigation options: "Back", "Alaris, Nancy Nh", and "Next". The main content area displays the medication name and concentration: "DOPamine 400 mg in dextrose 5 % 250 mL infusion (premix)", followed by "New Bag, 6/1/23 13:31". A green bar with white text says "Select a time and tap Next". Below this is a warning icon and text: "Order details will not be sent to the pump because a pump start time is selected." The medication ID "TLVP 15789008" is shown, followed by a "Channel Type" section with "Primary" selected and "Piggyback" as an option. The "Association Start" section has "Now" selected. The dose field shows "6/1/23 13:31 (3 mcg/kg/min)" selected, with "Other" as an alternative. The interface is clean with a light blue and white color scheme.

7. Click Next. The administration is now saved.

# Problem Solving

## **Basic Troubleshooting**

If you are not able to start or back associate an infusion with interop, first ask:

- Is the pump on?
- Is the order verified?
- Is the fluid/drug in scope?
  - Check to MAR for the “Manually Program Pump for this Medication” message
- Is the administration time correct?
- Is the drug ordered as an “infusion”?
  - Other IV routes (i.e. IV Push/Injection) are out of scope for interoperability.
- Does the order have a rate?
  - The order MUST have a rate. This can be seen on the MAR before sending order details.
- Was the infusion manually programmed in Basic Mode?
  - Basic infusions will need to be promoted to Guardrails before using interop. “Promoted” action can be found under OPTION keys.
- Is the pump already associated with a “running” order?
  - Each channel can only be associated with 1 infusion at a time. Nurse will be prompted to disassociate if the channel is already associated with another order.
- Is the pump ready?
  - The pump must be powered on, have the current drug library, in the appropriate profile, have WiFi, and have a barcode.




# Problem Solving (cont.)

## Error Messages

Warnings and Error messages will display with precise wording and direction. You should follow the directions on the error message and if it is an error that you can resolve, click “Resend Order”. Remember manual program pump means you are NOT using interoperability, so that should be selected only when you cannot resolve the issue.

(see page 8 #11)



- Error message: Medication ID Not Found
  - Technical issue that must be fixed by IT.
  - Move forward without interoperability and create a help desk ticket.
- Error message: Rate Required
  - Enter the rate in the MAR admin window and resend. Create a help desk ticket.

## Infusion Verify Troubleshooting

- Larger than expected volumes
  - This may happen if stop actions were not documented on the previous administration.
  - To fix: See page 22.
- Order associated, but not showing in the Infusion Verify window
  - Flowsheet rows must be active to document in Infusion Verify.

# Problem Solving (cont.)


- To fix: Check flowsheets to confirm that you can see the flowsheet rows for that infusion. If you cannot, select Unhide Completed to reveal, then complete Infusion Verify. Infusion groups must be checked in Flowsheets.

## **Alaris Pump Barcode:**

- If the barcode sticker falls off the pump channel (or will not scan), call Bio Med to get the sticker replaced and a new pump.
- Does the barcode match serial number? Options > Page Down>Serial Numbers

## **Association Error:**

If you receive a pump already associated message this is because the previous patient was not disassociated from the pump.



1. If the error is showing the incorrect patient, click Yes to continue with the disassociation.
2. Associate the pump with the patient to whom you will be administering the medications

## **Calculated Volume for an Offline Pump**

If the pump cannot communicate with Epic due to downtime, pump off, or a wireless issue, Infusion Verify will calculate the rate that was set before the interruption in communication.

# Problem Solving (cont.)

If the rate was manually changed during the interruption, Epic will show a higher or lower rate than what is showing on the pump to calculate volumes infused.

You will need to compare the rate on the pump with the rate showing in the Infusion Verify window in Epic, and make sure they match. If they do not match, you will have to manually enter the rate and volume infused using the physical pump as your guide.

## Time Lock State / Day Light Savings

If time change has occurred and you are unable to send order details to the pump, Epic is not communicating, or if information from an infusion is not appearing in Infusion Verify then the pump could be in a time lock state. Follow the steps below to check if it is locked.

With current infusions running:

1. Press Options on the Pump
2. Click the Page Down soft key.
3. Select Time of Day.
4. Validate that the time of day matches the expected time.
  - a. If it matches click EXIT.
  - b. If it does not match – Do NOT press change time or confirm.
    - POWER DOWN the entire system
    - Then press ON

\*This will re-sync the time AND correct the problem
  - c. Select No to new patient
  - d. Click Yes to confirm the same profile
  - e. Select the channel for each infusion
  - f. Click Restore to restore each infusion
  - g. Press Start

# Frequently Asked Questions

## **Can I clear my pump volumes?**

Yes, even though Epic is pulling all of your volumes, it is a good practice to routinely clear the volume infused on your pumps. This does not effect your I&O's in EPIC, but can be helpful – especially when infusing both “in” and “out” of scope medications.

## **I just received a patient from an “out of scope” department with medications infusing. The medications have all been charted in the MAR. What do I do?**

You would need to “Associate the pumps”. See page 19.

## **What happens during EPIC downtime?**

One-time medications that are started and stopped during downtime can not be back associated – these volumes would have to be charted manually in the MAR & I/O flowsheet.

If the medication was started during downtime AND still infusing, the medication can be back associated. (see page 20)

## **What if I temporarily loose WIFI (severe weather or transporting via an elevator)?**

Interoperability stores data for up to 8 hours. Once the pump gains WIFI, that information can be pulled over into EPIC via Infusion Verify.

# Frequently Asked Questions

**When I try to Infusion Verify, I do not see a medication or its associated volumes?**

There is a slight delay in transmitting data.

Example: Medication administered at 13:00:02 – data will not populate until the following “minute” (must wait 58 seconds)– 13:01:00

Medication administered at 13:00:58 – (must wait 2 seconds) –until the following “minute” – 13:01:00

If data still not pulling & WIFI is present, contact Help Desk to place a ticket.

**When I Infusion Verify, I notice some of my rows are “greyed-out”?**

This will appear when there are pauses that last less than 5 minutes.

**I only want to Verify one of my continuous infusions, can I do that?**

Yes, when you Infusion Verify, you can unclick any boxes that you do not wish to document. You can also click the medication in the I/O Flowsheet to verify that medication alone.

**I gave a medication emergently and I’m trying to “Back Associate”, but it is not asking me to “scan the pump”?**

Verify that there is WIFI and the order in in EPIC and verified by pharmacy.

If yes to above – a medication must be infusing for at least 5 minutes in order to “back associate”.

And do



# Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

# Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---